


Children and Young People's Diabetes Team

How to install Diasend at home


Diasend is now available for families to use to download information from insulin pumps, blood glucose meters and continuous glucose monitors at home. You will also be able to share this information with the diabetes team.

To download at home go to the Diasend website www.diasend.com




The screenshot shows the Diasend website homepage. At the top, there is a navigation menu with 'Clinic', 'Patient', 'Support', and 'About us'. Below the menu, a large banner features the text 'diasend now easily available for all patients' and a yellow circular button that says 'Click here to register for FREE!'. Below the banner, there are two video thumbnails: 'diasend film' and 'diasend Patient registration'. To the right of these thumbnails is a login form with fields for 'Username' and 'Password', a 'Log in' button, and a 'Forgot your password?' link. Below the login form is a yellow 'Register here' button. A green line points from the 'Register here' button to a callout box on the right that says 'Click here to go to patient registration'. At the bottom of the page, there is a 'Contact us' section and a footer with security information.

Either select the patient tab at the top of the page, or click on the orange circle, you will then see a screen that provides some information videos and allows you to set up an account to upload your devices. You will need to download the Diasend Uploader.




Click here to start the registration process



Choose the GB flag icon

You will then be able to set up a user name and password; you will need the Alder Hey clinic ID which is **37-65297** to allow your data to be linked to Alder Hey



The screenshot shows a web browser window displaying the registration page for diasend.com. The page is titled "Login details" and is part of a four-step process: 1. Login details, 2. Personal information, 3. Share data, and 4. Confirm registration. The main content area says "You're just a minute away from having your free diasend® Personal account" and "Login details". Below this, there are input fields for "E-mail" and "Repeat e-mail", "Password", and "Repeat password". A note states: "Your e-mail address will be your user name for diasend®, The e-mail address you enter must be valid." Another note says: "Your password must contain a minimum of 8 characters." The "diasend.com" logo is at the bottom left, and a "CONTINUE" button is at the bottom right. On the right side of the browser window, there is a sidebar with the "Hey" logo and system information: "AH015531", "10.17.56.41", "Windows 7", "3497 MB", "Service Pack 1", "8.0.7601.17514", "26/01/2015 08:53", and "Centre: 2343" with the URL "alderhey.nhs.uk".

If you need help setting up the Diasend System at home please contact a member of the diabetes team by email to diabetes@alderhey.nhs.uk or by calling your key worker.

Diasend Patient Registration – January 2015
Diabetes Team
Alder Hey
www.alderhey.nhs.uk
0151 228 4811